

Anti-Cd14 antibody [Sa14-2] {PE/Cyanine7} (STJA0000985)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	0.2 mg/mL
Conjugation	PE/Cyanine7
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2ak
Storage Instruction	Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	12475
Gene Symbol	Cd14
UniProt ID	CD14_MOUSE

ADDITIONAL INFORMATION

Note **FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.**



Mouse abdominal macrophages elicited by starch broth were stained with Anti-CD14 antibody (STJA0000985) [Used at 0.2 $\mu\text{g}/10^6$ cells dilution] (filled gray histogram). Unstained macrophages (blank black histogram) were used as control.