

## Anti-CD25/IL2RA antibody [BC96] {PE/Cyanine7} (STJA0001042)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | FC                 |
| Host / Source | Mouse              |
| Reactivity    | Human              |

### PRODUCT PROPERTIES


|                     |                                                                                                  |
|---------------------|--------------------------------------------------------------------------------------------------|
| Clonality           | Monoclonal                                                                                       |
| Concentration       | 5 Mu L/Test                                                                                      |
| Conjugation         | PE/Cyanine7                                                                                      |
| Formulation         | Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.      |
| Isotype             | IgG1k                                                                                            |
| Storage Instruction | Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product. |

### TARGET INFORMATION

|             |                             |
|-------------|-----------------------------|
| Gene ID     | <a href="#">3559</a>        |
| Gene Symbol | <a href="#">IL2RA</a>       |
| UniProt ID  | <a href="#">IL2RA_HUMAN</a> |

### ADDITIONAL INFORMATION

Note FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.



Human peripheral blood lymphocytes were stained with Anti-Human CD4 Monoclonal Antibody (APC Conjugated) and Anti-CD25 antibody (STJA0001042) (Left). Lymphocytes stained with Anti-Human CD4 Monoclonal Antibody (APC Conjugated) and Mouse IgG1 Isotype Control PE/Cy7 (Right) were used as control.