

Anti-CD11a/ITGAL antibody [R7-1] {PE/Cyanine5} (STJA0001050)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	5 Mu L/Test
Conjugation	PE/Cyanine5
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG1k
Storage Instruction	Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	3683
Gene Symbol	ITGAL
UniProt ID	ITAL_HUMAN

ADDITIONAL INFORMATION

Note **FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.**



Human peripheral blood lymphocytes were stained with Anti-CD11a antibody (STJA0001050) (filled gray histogram) or Mouse IgG1 Isotype Control PE/Cy5 (empty black histogram).