

## Anti-CD44H antibody [OX-49] {PE} (STJA0001070)

### GENERAL INFORMATION


Product Type	Primary antibodies
Applications	FC
Host / Source	Mouse
Reactivity	Rat

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	5 Mu L/Test
Conjugation	PE
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2ak
Storage Instruction	Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

### ADDITIONAL INFORMATION

Note **FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.**



Rat splenocytes were stained with Anti-CD44H antibody (STJA0001070) (filled gray histogram). Unstained splenocytes (empty black histogram) were used as control.