

Anti-CD44H antibody [OX-49] {PE/Cyanine5} (STJA0001077)

GENERAL INFORMATION


Product Type	Primary antibodies
Applications	FC
Host / Source	Mouse
Reactivity	Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	0.2 mg/mL
Conjugation	PE/Cyanine5
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2ak
Storage Instruction	Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

ADDITIONAL INFORMATION

Note FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.



Rat splenocytes were stained with Anti-CD44H antibody (STJA0001077) [Used at 0.2 μ g/ 10^6 cells dilution] (filled gray histogram). Unstained splenocytes (empty black histogram) were used as control.