

Anti-CD90.1 antibody [OX-7] {PE} (STJA0001083)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Mouse
Reactivity	Rat

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	0.2 mg/mL
Conjugation	PE
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG1k
Storage Instruction	Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	24832
Gene Symbol	Thy1
UniProt ID	THY1_RAT

ADDITIONAL INFORMATION

Note FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.



Rat thymocytes were stained with Anti-CD90/Mouse CD90.1 antibody (STJA0001083) [Used at 0.05 µg/10⁶ cells dilution] (filled gray histogram). Unstained thymocytes (empty black histogram) were used as control.