

## Anti-CD45/Ptprc antibody [OX-1] {APC} (STJA0001087)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Mouse
Reactivity	Rat

### PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	5 Mu L/Test
Conjugation	APC
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG1k
Storage Instruction	Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

### TARGET INFORMATION

Gene ID	<a href="#">24699</a>
Gene Symbol	<a href="#">Ptprc</a>
UniProt ID	<a href="#">PTPRC_RAT</a>

### ADDITIONAL INFORMATION

Note **FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.**



Rat splenocytes were stained with Anti-CD45 antibody (STJA0001087) (filled gray histogram). Unstained splenocytes (empty black histogram) were used as control.