

Anti-CD45/Ptprc antibody [OX-1] {PE/Cyanine5} (STJA0001094)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Mouse
Reactivity	Rat

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	0.2 mg/mL
Conjugation	PE/Cyanine5
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG1k
Storage Instruction	Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	24699
Gene Symbol	Ptprc
UniProt ID	PTPRC_RAT

ADDITIONAL INFORMATION

Note FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.



Rat splenocytes were stained with Anti-CD45 antibody (STJA0001094) [Used at 0.2 µg/10⁶ cells dilution] (filled gray histogram). Unstained splenocytes (empty black histogram) were used as control.