

## Anti-CD3 antibody [G4.18] {FITC} (STJA0001101)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Mouse
Reactivity	Rat

### PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	0.5 mg/mL
Conjugation	FITC
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG3k
Storage Instruction	Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

### TARGET INFORMATION

Gene ID	<a href="#">25710</a> <a href="#">300678</a>
Gene Symbol	<a href="#">Cd3d</a> <a href="#">Cd3g</a>
UniProt ID	<a href="#">CD3D_RAT</a> <a href="#">CD3G_RAT</a>

### ADDITIONAL INFORMATION

Note FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.



Rat splenocytes were stained with Anti-CD3 antibody (STJA0001101) [Used at 0.2  $\mu$ g/ $10^6$  cells dilution] (filled gray histogram). Unstained splenocytes (empty black histogram) were used as control.