

Anti-Foxp3 antibody [3G3] {APC} (STJA0001118)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ICFC
Host / Source	Mouse
Reactivity	Mouse

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	5 Mu L/Test
Conjugation	APC
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG1k
Storage Instruction	Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	50943 20371
Gene Symbol	FOXP3 Foxp3
UniProt ID	FOXP3_HUMAN FOXP3_MOUSE

ADDITIONAL INFORMATION

Note FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.



C57BL/6 splenocytes were surface stained with Anti-Mouse CD4 Monoclonal Antibody (FITC Conjugated) followed by fixation and permeabilization using the Foxp3 Staining/Transcription Factor Buffer Set and intracellular staining with Anti-Foxp3 antibody (STJA0001118).