

Anti-MHC 2 (I-A/I-E) antibody [M5/114] (STJA0001218)

GENERAL INFORMATION	
Product Type	Primary antibodies
Applications	Block/FA/FC/WB/IP/IF
Host / Source	Rat
Reactivity	Mouse
PRODUCT PROPERTIES	
Clonality	Monoclonal
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Formulation	0.2 μ m filtered in PBS, pH7.2. Azide Free (AF) /Low Endotoxin (LE) : Contains no stabilizers or preservatives. Endotoxin level is < 2 EU/mg as determined by LAL gel clotting assay.
Isotype	IgG2bk
Storage Instruction	This product can be stored at 2-8°C for 12 months. Do not freeze.
TARGET INFORMATION	
Gene ID	14961
Gene Symbol	H2-Ab1
UniProt ID	HB2A_MOUSE
ADDITIONAL INFORMATION	
Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.