

## Anti-Cd80 antibody [16-10A1] (STJA0001228)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC/IHC/IP
Host / Source	Armenian Hamster
Reactivity	Mouse

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Formulation	PBS with 0.05% Proclin300
Isotype	IgG
Storage Instruction	This product can be stored at 2-8°C for 12 months. Do not freeze.

### TARGET INFORMATION

Gene ID	<a href="#">12519</a>
Gene Symbol	<a href="#">Cd80</a>
UniProt ID	<a href="#">CD80_MOUSE</a>

### ADDITIONAL INFORMATION

Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.
------	---