

Anti-Cd86 antibody [GL-1] (STJA0001234)

GENERAL INFORMATION

| | |
|---------------|--------------------|
| Product Type | Primary antibodies |
| Applications | FC/IHC-F/IP/ICC |
| Host / Source | Rat |
| Reactivity | Mouse |

PRODUCT PROPERTIES

| | |
|---------------------|---|
| Clonality | Monoclonal |
| Concentration | 0.5 mg/mL |
| Conjugation | Unconjugated |
| Formulation | 0.2 μ m filtered in PBS, pH7.2. Azide Free (AF) /Low Endotoxin (LE) : Contains no stabilizers or preservatives. Endotoxin level is < 2 EU/mg as determined by LAL gel clotting assay. |
| Isotype | IgG2ak |
| Storage Instruction | This product can be stored at 2-8°C for 12 months. Do not freeze. |

TARGET INFORMATION

| | |
|-------------|----------------------------|
| Gene ID | 12524 |
| Gene Symbol | Cd86 |
| UniProt ID | CD86_MOUSE |

ADDITIONAL INFORMATION

| | |
|------|---|
| Note | FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS. |
|------|---|