

Anti-CD16/32 antibody [2.4G2] {AF488} (STJA0001247)

GENERAL INFORMATION

| | |
|---------------|--------------------|
| Product Type | Primary antibodies |
| Applications | FC |
| Host / Source | Rat |
| Reactivity | Mouse |

PRODUCT PROPERTIES

| | |
|---------------------|--|
| Clonality | Monoclonal |
| Concentration | 5 Mu L/Test |
| Conjugation | AF488 |
| Formulation | Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant. |
| Isotype | IgG2bk |
| Storage Instruction | Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product. |

TARGET INFORMATION

| | |
|-------------|--|
| Gene ID | 14130 NA |
| Gene Symbol | Fcgr2 Fcgr3 |
| UniProt ID | FCGR2_MOUSE FCGR3_MOUSE |

ADDITIONAL INFORMATION

| | |
|------|---|
| Note | FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS. |
|------|---|