

**Anti-CD62L/Sell antibody [Me114] {PE/Cyanine5.5} (STJA0001269)****GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	FC
Host / Source	Rat
Reactivity	Mouse

**PRODUCT PROPERTIES**

Clonality	Monoclonal
Concentration	0.2 mg/mL
Conjugation	PE/Cyanine5.5
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2ak
Storage Instruction	Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

**TARGET INFORMATION**

Gene ID	<a href="#">20343</a>
Gene Symbol	<a href="#">Sell</a>
UniProt ID	<a href="#">LYAM1_MOUSE</a>

**ADDITIONAL INFORMATION**

Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.
------	---