

Anti-CD3 antibody [17A2] (STJA0001277)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC/IHC/IP
Host / Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Formulation	PBS with 0.05% Proclin300
Isotype	IgG2bk
Storage Instruction	This product can be stored at 2-8°C for 12 months. Do not freeze.

TARGET INFORMATION

Gene ID	12501 12502 12500 Cd3e Cd3g CD3E_MOUSE CD3G_MOUSE CD3D_MOUSE
---------	---

ADDITIONAL INFORMATION

Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.
------	---