

Anti-CD3 antibody [17A2] {Biotin} (STJA0001278)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	0.5 mg/mL
Conjugation	Biotin
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2bk
Storage Instruction	This product can be stored at 2-8°C for 12 months. Do not freeze.

TARGET INFORMATION

Gene ID	12500 12501 12503 Cd3d Cd3e CD3D_MOUSE CD3E_MOUSE CD3Z_MOUSE
---------	---

ADDITIONAL INFORMATION

Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.
------	---