

## Anti-CD132/Il2rg antibody [3E12] {Biotin} (STJA0001283)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Rat
Reactivity	Mouse

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	0.5 mg/mL
Conjugation	Biotin
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2bk
Storage Instruction	This product can be stored at 2-8°C for 12 months. Do not freeze.

### TARGET INFORMATION

Gene ID	<a href="#">16186</a>
Gene Symbol	<a href="#">Il2rg</a>
UniProt ID	<a href="#">IL2RG_MOUSE</a>

### ADDITIONAL INFORMATION

Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.
------	-------------------------------------------------------------------------------------------------------------------