

## Anti-CD47 antibody [MIAP410] (STJA0001289)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Mouse
Reactivity	Human/Mouse/Rat

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Formulation	PBS with 0.05% Proclin300
Isotype	IgG1k
Storage Instruction	This product can be stored at 2-8°C for 12 months. Do not freeze.

### TARGET INFORMATION

Gene ID	<a href="#">961</a> <a href="#">29364</a> <a href="#">16423</a> <a href="#">CD47</a> <a href="#">Cd47</a> <a href="#">CD47_HUMAN</a> <a href="#">CD47_RAT</a> <a href="#">CD47_MOUSE</a>
---------	---

### ADDITIONAL INFORMATION

Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.
------	---