

Anti-Cd40 antibody [FGK4.5/FGK45] {PE/Cyanine7} (STJA0001356)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	FC
Host / Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	5 Mu L/Test
Conjugation	PE/Cyanine7
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2ak
Storage Instruction	Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	21939
Gene Symbol	Cd40
UniProt ID	TNR5_MOUSE

ADDITIONAL INFORMATION

Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.
------	---