

## Anti-CD162/Selplg antibody [4RA10] (STJA0001397)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | FC/WB/Block/IP     |
| Host / Source | Rat                |
| Reactivity    | Mouse              |

### PRODUCT PROPERTIES

|                     |                                                                   |
|---------------------|-------------------------------------------------------------------|
| Clonality           | Monoclonal                                                        |
| Concentration       | 0.5 mg/mL                                                         |
| Conjugation         | Unconjugated                                                      |
| Formulation         | PBS with 0.05% Proclin300                                         |
| Isotype             | IgG1k                                                             |
| Storage Instruction | This product can be stored at 2-8°C for 12 months. Do not freeze. |

### TARGET INFORMATION

|             |                                                                  |
|-------------|------------------------------------------------------------------|
| Gene Symbol | <a href="#">Selplg</a> Uniprot ID= <a href="#">SELPL_MOUSE</a> " |
|-------------|------------------------------------------------------------------|

### ADDITIONAL INFORMATION

|      |                                                                                                                   |
|------|-------------------------------------------------------------------------------------------------------------------|
| Note | FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------|