

Anti-CD120b/Tnfrsf1b antibody [TR75-54.7] (STJA0001414)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC/WB/IP
Host / Source	Armenian Hamster
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Formulation	PBS with 0.05% Proclin300
Isotype	IgG
Storage Instruction	This product can be stored at 2-8°C for 12 months. Do not freeze.

TARGET INFORMATION

Gene ID	21938
Gene Symbol	Tnfrsf1b
UniProt ID	TNR1B_MOUSE

ADDITIONAL INFORMATION

Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.
------	---