

**Anti-CD120b/Tnfrsf1b antibody [TR75-54.7] {Biotin} (STJA0001415)****GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	FC
Host / Source	Armenian Hamster
Reactivity	Mouse

**PRODUCT PROPERTIES**

Clonality	Monoclonal
Concentration	0.5 mg/mL
Conjugation	Biotin
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG
Storage Instruction	Store at 2-8°C. Do not freeze this product.

**TARGET INFORMATION**

Gene ID	<a href="#">21938</a>
Gene Symbol	<a href="#">Tnfrsf1b</a>
UniProt ID	<a href="#">TNR1B_MOUSE</a>

**ADDITIONAL INFORMATION**

Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.
------	---