

## Anti-CD57/B3GAT1 antibody [HNK-1] (STJA0001485)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC/WB/IHC
Host / Source	Mouse
Reactivity	Human

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Formulation	PBS with 0.05% Proclin300
Isotype	IgMk
Storage Instruction	This product can be stored at 2-8°C for 12 months. Do not freeze.

### TARGET INFORMATION

Gene ID	<a href="#">27087</a>
Gene Symbol	<a href="#">B3GAT1</a>
UniProt ID	<a href="#">B3GA1_HUMAN</a>

### ADDITIONAL INFORMATION

Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.
------	---