

## Anti-CD15/FUT4 antibody [HI98] (STJA0001514)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC/IHC
Host / Source	Mouse
Reactivity	Human

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Formulation	PBS with 0.05% Proclin300
Isotype	IgMk
Storage Instruction	This product can be stored at 2-8°C for 12 months. Do not freeze.

### ADDITIONAL INFORMATION

Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.
------	---