

Anti-CD11b antibody [M1/70] (STJA0001524)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC/IHC/IP
Host / Source	Rat
Reactivity	Mouse/Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Formulation	PBS with 0.05% Proclin300
Isotype	IgG2bk
Storage Instruction	This product can be stored at 2-8°C for 12 months. Do not freeze.

TARGET INFORMATION

Gene ID	3684 NA
Gene Symbol	ITGAM Itgam
UniProt ID	ITAM_HUMAN ITAM_MOUSE

ADDITIONAL INFORMATION

Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.
------	---