

**Anti-CD103/Itgae antibody [M290] {PE/Cyanine5.5} (STJA0001557)****GENERAL INFORMATION**

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | FC                 |
| Host / Source | Rat                |
| Reactivity    | Mouse              |

**PRODUCT PROPERTIES**

|                     |                                                                                                  |
|---------------------|--------------------------------------------------------------------------------------------------|
| Clonality           | Monoclonal                                                                                       |
| Concentration       | 5 Mu L/Test                                                                                      |
| Conjugation         | PE/Cyanine5.5                                                                                    |
| Formulation         | Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.      |
| Isotype             | IgG2ak                                                                                           |
| Storage Instruction | Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product. |

**TARGET INFORMATION**

|             |                            |
|-------------|----------------------------|
| Gene ID     | <a href="#">16407</a>      |
| Gene Symbol | <a href="#">Itgae</a>      |
| UniProt ID  | <a href="#">ITAE_MOUSE</a> |

**ADDITIONAL INFORMATION**

|      |                                                                                                                   |
|------|-------------------------------------------------------------------------------------------------------------------|
| Note | FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------|