

## Anti-CD90.1 antibody [19E12] (STJA0001602)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | FC                 |
| Host / Source | Rat                |
| Reactivity    | Mouse              |

### PRODUCT PROPERTIES

|                     |   |
|---------------------|---|
| Clonality           | Monoclonal  |
| Concentration       | 0.5 mg/mL   |
| Conjugation         | Unconjugated  |
| Formulation         | PBS with 0.05% Proclin300   |
| Isotype             | IgG2ak  |
| Storage Instruction | This product can be stored at 2-8°C for 12 months. Do not freeze. |

### TARGET INFORMATION

|             |                            |
|-------------|----------------------------|
| Gene ID     | <a href="#">21838</a>      |
| Gene Symbol | <a href="#">Thy1</a>       |
| UniProt ID  | <a href="#">THY1_MOUSE</a> |

### ADDITIONAL INFORMATION

|      |   |
|------|---|
| Note | FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS. |
|------|---|