

Anti-Cd8a antibody [OX-8] {PE/Cyanine7} (STJA0001626)

GENERAL INFORMATION		
	Product Type	Primary antibodies
	Applications	FC
	Host / Source	Mouse
	Reactivity	Rat
PRODUCT PROPERTIES		
	Clonality	Monoclonal
	Concentration	5 Mu L/Test
	Conjugation	PE/Cyanine7
	Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
	Isotype	IgG1k
	Storage Instruction	Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.
TARGET INFORMATION		
	Gene ID	24930
	Gene Symbol	Cd8a
	UniProt ID	CD8A_RAT
ADDITIONAL INFORMATION		
	Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.