

Anti-Ly6g antibody [1A8] (STJA0001695)

GENERAL INFORMATION		
	Product Type	Primary antibodies
	Applications	Depletion/FC/IHC/IF
	Host / Source	Rat
	Reactivity	Mouse
PRODUCT PROPERTIES		
	Clonality	Monoclonal
	Concentration	0.5 mg/mL
	Conjugation	Unconjugated
	Formulation	0.2 µm filtered in PBS, pH7.2. Azide Free (AF) /Low Endotoxin (LE) : Contains no stabilizers or preservatives. Endotoxin level is < 2 EU/mg as determined by LAL gel clotting assay.
	Isotype	IgG2ak
	Storage Instruction	This product can be stored at 2-8°C for 12 months. Do not freeze.
TARGET INFORMATION		
	Gene ID	546644
	Gene Symbol	Ly6g
	UniProt ID	LY6G_MOUSE
ADDITIONAL INFORMATION		
	Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.