

**Anti-CD45R/B220 antibody [RA3.3A 1/6.1] (STJA0001711)****GENERAL INFORMATION**

|               |                                  |
|---------------|----------------------------------|
| Product Type  | Primary antibodies               |
| Applications  | Cell Sep-Neg/Depletion/FC/IHC/IP |
| Host / Source | Rat                              |
| Reactivity    | Mouse                            |

**PRODUCT PROPERTIES**

|                     |                                                                                                                                                                                           |
|---------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Clonality           | Monoclonal                                                                                                                                                                                |
| Concentration       | 0.5 mg/mL                                                                                                                                                                                 |
| Conjugation         | Unconjugated                                                                                                                                                                              |
| Formulation         | 0.2 $\mu$ m filtered in PBS, pH7.2. Azide Free (AF) /Low Endotoxin (LE) : Contains no stabilizers or preservatives. Endotoxin level is < 2 EU/mg as determined by LAL gel clotting assay. |
| Isotype             | IgMk                                                                                                                                                                                      |
| Storage Instruction | This product can be stored at 2-8°C for 12 months. Do not freeze.                                                                                                                         |

**TARGET INFORMATION**

|             |                                                            |
|-------------|------------------------------------------------------------|
| Gene ID     | <a href="#">19264</a><br><a href="#">5788</a>              |
| Gene Symbol | <a href="#">Ptprc</a><br><a href="#">PTPRC</a>             |
| UniProt ID  | <a href="#">PTPRC_MOUSE</a><br><a href="#">PTPRC_HUMAN</a> |

**ADDITIONAL INFORMATION**

|      |                                                                                                                   |
|------|-------------------------------------------------------------------------------------------------------------------|
| Note | FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------|