

Anti-CD45R/B220 antibody [RA3.3A 1/6.1] (STJA0001712)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC/IHC/IP
Host / Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Formulation	PBS with 0.05% Proclin300
Isotype	IgMk
Storage Instruction	This product can be stored at 2-8°C for 12 months. Do not freeze.

TARGET INFORMATION

Gene ID	5788 19264
Gene Symbol	PTPRC Ptprc
UniProt ID	PTPRC_HUMAN PTPRC_MOUSE

ADDITIONAL INFORMATION

Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.
------	---