

Anti-Ly6C antibody [Monts 1] (STJA0001758)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	Depletion/FC
Host / Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Formulation	0.2 μ m filtered in PBS, pH7.2. Azide Free (AF) /Low Endotoxin (LE) : Contains no stabilizers or preservatives. Endotoxin level is < 2 EU/mg as determined by LAL gel clotting assay.
Isotype	IgG2ak
Storage Instruction	This product can be stored at 2-8°C for 12 months. Do not freeze.

TARGET INFORMATION

Gene ID	17067 100041546
Gene Symbol	Ly6c1 Ly6c2
UniProt ID	LY6C1_MOUSE LY6C2_MOUSE

ADDITIONAL INFORMATION

Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.
------	---