

**Anti-CD45.2/Ptprc antibody [104.2] (STJA0001761)****GENERAL INFORMATION**

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | Block/FC/IHC       |
| Host / Source | Mouse              |
| Reactivity    | Mouse              |

**PRODUCT PROPERTIES**

|                     |                                                                                                                                                                                           |
|---------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Clonality           | Monoclonal                                                                                                                                                                                |
| Concentration       | 0.5 mg/mL                                                                                                                                                                                 |
| Conjugation         | Unconjugated                                                                                                                                                                              |
| Formulation         | 0.2 $\mu$ m filtered in PBS, pH7.2. Azide Free (AF) /Low Endotoxin (LE) : Contains no stabilizers or preservatives. Endotoxin level is < 2 EU/mg as determined by LAL gel clotting assay. |
| Isotype             | IgG2ak                                                                                                                                                                                    |
| Storage Instruction | This product can be stored at 2-8°C for 12 months. Do not freeze.                                                                                                                         |

**TARGET INFORMATION**

|             |                             |
|-------------|-----------------------------|
| Gene ID     | <a href="#">19264</a>       |
| Gene Symbol | <a href="#">Ptprc</a>       |
| UniProt ID  | <a href="#">PTPRC_MOUSE</a> |

**ADDITIONAL INFORMATION**

|      |                                                                                                                   |
|------|-------------------------------------------------------------------------------------------------------------------|
| Note | FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------|