

## Anti-CD19 antibody [4G7] {PE/TR} (STJA0001805)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | FC                 |
| Host / Source | Mouse              |
| Reactivity    | Human              |

### PRODUCT PROPERTIES

|                     |                                                                                             |
|---------------------|---------------------------------------------------------------------------------------------|
| Clonality           | Monoclonal                                                                                  |
| Concentration       | 5 Mu L/Test                                                                                 |
| Conjugation         | PE/TR                                                                                       |
| Purification        | Affinity purified                                                                           |
| Formulation         | Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant. |
| Isotype             | IgG1k                                                                                       |
| Storage Instruction | This product can be stored at 2-8°C for 12 months. Do not freeze.                           |

### TARGET INFORMATION

|             |                            |
|-------------|----------------------------|
| Gene ID     | <a href="#">930</a>        |
| Gene Symbol | <a href="#">CD19</a>       |
| UniProt ID  | <a href="#">CD19_HUMAN</a> |

### ADDITIONAL INFORMATION

|      |                                                                                                                   |
|------|-------------------------------------------------------------------------------------------------------------------|
| Note | FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------|