

Anti-CD279/PDCD1 antibody [J116] (STJA0001810)

GENERAL INFORMATION

| | |
|---------------|--------------------|
| Product Type | Primary antibodies |
| Applications | FC/IHC/IP |
| Host / Source | Mouse |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|---------------------|---|
| Clonality | Monoclonal |
| Concentration | 0.2 mg/mL |
| Conjugation | Unconjugated |
| Formulation | PBS with 0.05% Proclin300 |
| Isotype | IgG1k |
| Storage Instruction | This product can be stored at 2-8°C for 12 months. Do not freeze. |

TARGET INFORMATION

| | |
|-------------|-----------------------------|
| Gene ID | 5133 |
| Gene Symbol | PDCD1 |
| UniProt ID | PDCD1_HUMAN |

ADDITIONAL INFORMATION

| | |
|------|---|
| Note | FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS. |
|------|---|