

Anti-Cd274/PD-L1 antibody [10F.9G2] (STJA0001821)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	FC/IHC/IF
Host / Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Formulation	PBS with 0.05% Proclin300
Isotype	IgG2bk
Storage Instruction	This product can be stored at 2-8°C for 12 months. Do not freeze.

TARGET INFORMATION

Gene ID	60533
Gene Symbol	Cd274
UniProt ID	PD1L1_MOUSE

ADDITIONAL INFORMATION

Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.
------	---