

Anti-Ly6a antibody [D7] {AF647} (STJA0002086)

GENERAL INFORMATION		
	Product Type	Primary antibodies
	Applications	FC
	Host / Source	Rat
	Reactivity	Mouse
PRODUCT PROPERTIES		
	Clonality	Monoclonal
	Concentration	5 Mu L/Test
	Conjugation	AF647
	Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
	Isotype	IgG2ak
	Storage Instruction	Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.
TARGET INFORMATION		
	Gene ID	110454
	Gene Symbol	Ly6a
	UniProt ID	LY6A_MOUSE
ADDITIONAL INFORMATION		
	Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.