

## Anti-Il17a antibody [TC11-18H10.1] {PE} (STJA0002137)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ICFC
Host / Source	Rat
Reactivity	Mouse

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	0.2 mg/mL
Conjugation	PE
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG1k
Storage Instruction	Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

### TARGET INFORMATION

Gene ID	<a href="#">16171</a>
Gene Symbol	<a href="#">Il17a</a>
UniProt ID	<a href="#">IL17_MOUSE</a>

### ADDITIONAL INFORMATION

Note **FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.**

