

Anti-IL21 antibody [3A3-N2] {PE} (STJA0002150)

GENERAL INFORMATION

| | |
|---------------|--------------------|
| Product Type | Primary antibodies |
| Applications | ICFC |
| Host / Source | Mouse |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|---------------------|--|
| Clonality | Monoclonal |
| Concentration | 5 Mu L/Test |
| Conjugation | PE |
| Formulation | Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant. |
| Isotype | IgG1k |
| Storage Instruction | Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product. |

TARGET INFORMATION

| | |
|-------------|----------------------------|
| Gene ID | 59067 |
| Gene Symbol | IL21 |
| UniProt ID | IL21_HUMAN |

ADDITIONAL INFORMATION

Note FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.

