

## Anti-IL6 antibody [MP5-20F3] {APC} (STJA0002177)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ICFC
Host / Source	Rat
Reactivity	Mouse

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	0.2 mg/mL
Conjugation	APC
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG1k
Storage Instruction	Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

### TARGET INFORMATION

Gene ID	<a href="#">16193</a>
Gene Symbol	<a href="#">IL6</a>
UniProt ID	<a href="#">IL6_MOUSE</a>

### ADDITIONAL INFORMATION

Note FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.

