

Anti-CD44 antibody [Hermes-1] (STJA0002208)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC/IF
Host / Source	Rat
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Formulation	PBS with 0.05% Proclin300
Isotype	IgG2ak
Storage Instruction	This product can be stored at 2-8°C for 12 months. Do not freeze.

TARGET INFORMATION

Gene ID	960
Gene Symbol	CD44
UniProt ID	CD44_HUMAN

ADDITIONAL INFORMATION

Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.
------	---