

Anti-Foxp3 antibody [3G3] {FITC} (STJA0002249)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ICFC
Host / Source	Mouse
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	5 Mu L/Test
Conjugation	FITC
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG1k
Storage Instruction	Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	20371 50943
Gene Symbol	Foxp3 FOXP3
UniProt ID	FOXP3_MOUSE FOXP3_HUMAN

ADDITIONAL INFORMATION

Note FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.

