

Syrian Hamster IgG Isotype Control [SHG-1] {Biotin} (STJA0002272)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Syrian Hamster

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	0.5 mg/mL
Conjugation	Biotin
Formulation	PBS with 0.05% Proclin300, 1% BSA
Isotype	IgG
Storage Instruction	This product can be stored at 2-8°C for 12 months. Do not freeze.

ADDITIONAL INFORMATION

Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.
------	---