

## Syrian Hamster IgG Isotype Control [SHG-1] {AF488} (STJA0002282)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | FC                 |
| Host / Source | Syrian Hamster     |

### PRODUCT PROPERTIES

|                     |                                                                                                  |
|---------------------|--------------------------------------------------------------------------------------------------|
| Clonality           | Monoclonal                                                                                       |
| Concentration       | 0.5 mg/mL                                                                                        |
| Conjugation         | AF488                                                                                            |
| Formulation         | PBS with 0.05% Proclin300, 1% BSA                                                                |
| Isotype             | IgG                                                                                              |
| Storage Instruction | Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product. |

### ADDITIONAL INFORMATION

|      |                                                                                                                   |
|------|-------------------------------------------------------------------------------------------------------------------|
| Note | FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------|