

## Rat IgG1k Isotype Control [HRPN] {PE/Cyanine7} (STJA0002414)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Rat

### PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	PE/Cyanine7
Formulation	PBS with 0.05% Proclin300, 1% BSA
Isotype	IgG1k
Storage Instruction	Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

### ADDITIONAL INFORMATION

Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.
------	---