

**Anti-CD3/CD16/CD8a antibody [OKT3, 3G8, 5.1H11] {FITC/PE} (STJA0002509)****GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	FC
Reactivity	Human

**PRODUCT PROPERTIES**

Clonality	Monoclonal
Concentration	5 Mu L/Test
Conjugation	FITC/PE
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Storage Instruction	Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

**ADDITIONAL INFORMATION**

Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.
------	---