

**Anti-CD3/CD19 antibody [OKT3, 4G7] {FITC/PE} (STJA0002512)****GENERAL INFORMATION**

|              |                    |
|--------------|--------------------|
| Product Type | Primary antibodies |
| Applications | FC                 |
| Reactivity   | Human              |

**PRODUCT PROPERTIES**

|                     |                                                                                                  |
|---------------------|--------------------------------------------------------------------------------------------------|
| Clonality           | Monoclonal                                                                                       |
| Concentration       | 5 Mu L/Test                                                                                      |
| Conjugation         | FITC/PE                                                                                          |
| Formulation         | Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.      |
| Storage Instruction | Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product. |

**ADDITIONAL INFORMATION**

|      |                                                                                                                   |
|------|-------------------------------------------------------------------------------------------------------------------|
| Note | FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------|