

Anti-CD18/Itgb2 antibody [M18/2] {PE/Cyanine7} (STJA0002569)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	5 Mu L/Test
Conjugation	PE/Cyanine7
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2ak
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	16414
Gene Symbol	Itgb2
UniProt ID	ITB2_MOUSE

ADDITIONAL INFORMATION

Note FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.

