

Mouse IgG3k Isotype Control [A112-3] {AF488} (STJA0002747)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	0.5 mg/mL
Conjugation	AF488
Formulation	PBS with 0.05% Proclin300, 1% BSA
Isotype	IgG3k
Storage Instruction	Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

ADDITIONAL INFORMATION

Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.
------	---